

Tubing Cutters



Rotary Tubing Cutter

- Cuts 1/32" to 1/8" o.d. tubing burr-free
- Leaves tubing i.d. open
- Works with metal and glass-lined tubing



Heavy-duty Tubing Cutter

- Inexpensive, well-built tool
- Cuts 1/8" to 1-1/8" o.d. tubing
- Works with metal tubing




Plier-type Tubing Cutters

- Cuts 1/16" o.d. tubing
- Ideal for cutting tubing in tight places
- Cuts metal or plastic tubing



Ceramic Fused Silica Cutter

- Cleanly cuts fused silica tubing
- Small enough to keep handy in your pocket
- Just score and snap



Electric Tubing Cutter

- Cut 1/16", 1/8", and 1/4" o.d. tubings with internal diameters as small as 0.008"
- Requires no lubricant or cutting fluid
- Assures square and burr-free cuts



Clean-cut Tubing Cutter

- Use with PTFE, Tefzel®, Polyolefin, PEEK
- Make clean, square cuts without deforming inner diameters
- Cuts 1/32" to 1/8" o.d.



Plastic Tubing Cutter

- Ideal for PEEK, Tefzel®, and other plastic tubing
- Cuts any size up to 1/4" o.d.



Tubing Accessories

Flex-wrap Tubing

- Holds parallel tubing or wiring together
- Ideal for gradient pumps or purging-gas lines
- Chemically inert 100% PTFE



Flex-Wrap Tubing

Color	1/8"–1/2" Part No.
Blue	35913
Orange	35914
Black	35916
White	35917

Tubing Cutters

Description	Qty.	Part No.
Rotary Tubing Cutter	ea	3165
Replacement Cutting Wheels	3	3167
Heavy-Duty Tubing Cutter	ea	13995
Plier-Type Tubing Cutter	ea	3169
Ceramic Fused Silica Cutter	5	3194
Electric Tubing Cutter 110–120/220–240VAC, 50/60Hz	ea	10-0200
Replacement Parts and Accessories		
Dressing Tool for 1/8" and 1/4" o.d. Tubing	ea	10-0118
Replacement Cutting Wheels	3	10-0124
Replacement Dressing Tool for 1/16" o.d. Tubing	ea	10-0116
Clean-Cut™ Tubing Cutter	ea	35902
Replacement Blade	ea	35903
Plastic Tubing Cutter	ea	3206
Replacement Blade	ea	3214

tubing

Organizers

HPLC Columns and Accessories Organizers

- Holds 10 HPLC columns
- 17 compartments
- PVC construction
- 12" H x 24" W x 6" D
(30.48cm H x 60.96cm W x 15.24cm D)



HPLC Columns Organizers

- 30 HPLC columns
- Wall-mounts with two screws
- PVC construction with acrylic doors
- 17-3/8" H x 15" W x 2-7/8" D
(44.18cm H x 38.1cm W x 7.3cm D)



best seller

Plastic Box with Removable Dividers

- 13 removable dividers
- Ideal for fittings and bore-cut tubing



Lab Organizer with Adjustable Bins

- 14 adjustable dividers, 4 shelves
- PVC construction
- 20" H x 24" W x 10.5" D
(50.8cm H x 60.9cm W x 26.7cm D)



Organizer Box

- 9-compartment organizer includes eight 1-3/8" x 1-5/8" (3.5cm x 4.13cm) slots and one 1-1/8" x 3-1/4" (2.8cm x 8.2cm) slot
- 18-compartment organizer includes twelve 2" x 2" (5cm x 5cm) slots and six 2" x 4-3/16" (5cm x 10.6cm) slots



Organizers

Description	Part No.
HPLC Supplies Organizer	5238
HPLC Column Storage Unit	5240
Plastic Storage Box with Removable Dividers	803880
Lab Organizer with Adjustable Bins	5242
Organizer Box	
9-Compartment Box	6264
18-Compartment Box	6266

related products

Looking for vial organizers and storage racks?
See pages 380 and 381.



Tools

Magnifier

- For inspecting capillary column ends



20X Magnifier

Capillary Tool Kit

- Kit includes diamond scribe, two open end wrenches, PTFE-coated forceps, ruler, ferrule pick, 10' high temperature string, 20X magnifier, and pipe cleaner



5179

23-Piece Tool Kit in Zippered Case

- Kit includes seven screwdrivers, adjustable wrench, two pliers, wire stripper, knife, alignment tool, stainless rule, hex key set, scissors, two flexible files, burnisher, miniature soldering iron, solder aid, coil of solder, and desoldering braid



5323

Adjustable Wrench

- Chrome-Plated



5322

Open-End Wrench

- Chrome-Plated



5321

Handy-Lok Fitting Wrench

- Ideal for connections on injection valves



5641

9-Way Wrench

- For use with gas cylinders



5178

Miniature Files, Fine Cut

- Set includes six different-shaped, 5-1/2" long files
- Made of hard chrome steel suitable for smoothing and deburring



6279

Miniature Reamers

- Set includes six different sizes: 0.028"/0.018", 0.040"/0.026", 0.048"/0.032", 0.080"/0.050", 0.085"/0.053", 0.104"/0.062"



5324

Flame Checker

- Pivoting mirror head fogs when flame is lit



5715

Tools

Description	Qty.	Part No.
Magnifier, 20X	ea	3163
Capillary Tool Kit	ea	16879
23-Piece Tool Kit and Case	ea	16533
Adjustable Wrench		
4" Long—1/2" Capacity	ea	16549
6" Long—3/4" Capacity	ea	3211
8" Long—7/8" Capacity	ea	16522
Open-End Wrench		
1/4" x 5/16"—Fits 1/16" Fitting Nuts	ea	1998
7/16" x 1/2"—Fits 1/8" and 3/16" Fitting Nuts	ea	1999
9/16" x 5/8"—Fits 1/4" and 3/8" Fitting Nuts	ea	2000
Handy-Lok Fitting Wrench	ea	35751
9-Way Cylinder Wrench	ea	412933
Miniature Files, Fine Cut	6/pk	7237
Miniature Reamers	6/pk	7330
Flame Checker	ea	7047