


## Alltech® Solvent Filtration Apparatus

The Alltech® Solvent Filtration Apparatus rapidly filters and degasses HPLC and IC mobile phases. Depending on your analytical requirements, several interchangeable funnel and flask options are available. You can easily replace the funnel with the solvent pickup adapter. The Solvent Pickup Adapter draws solvent directly from a reagent bottle through a 47mm filter and into the collecting flask.



6266



2771

### Alltech® Solvent Filtration Apparatus

Description	Part No.
<i>Solvent Filtration Apparatus</i>	
1L Flask, 300mL Funnel, 47mm Filter Support, 47mm Clamp, 0.2µm Nylon Filters 100/pk	<b>2001</b>
2L Flask, 500mL Funnel, 47mm Filter Support, 47mm Clamp, 0.2µm Nylon Filters 100/pk	<b>2003</b>
<i>Replacement Parts and Accessories</i>	
1L Flask (For 2001)	<b>350062</b>
2L Flask (For 2003)	<b>350066</b>
300mL Glass Funnel (For 2001)	<b>350064</b>
500mL Glass Funnel (For 2003)	<b>350067</b>
1L Glass Funnel (for 2003)	<b>350070</b>
47mm Fritted Base Support (For 2001)	<b>350065</b>
47mm Fritted Base Support (For 2003)	<b>350068</b>
47mm Aluminum Clamp (For 2001, 2003)	<b>350061</b>
47mm Nylon Filter Membranes, 0.2µm, 100/pk	<b>2034</b>